


# BETACHEK *Visual*

## Blood Glucose Test Strips



BETACHEK Visual is a pink/blue “visually read” blood glucose test. BETACHEK Visual is a premium visual test with twin reactive pads and a 28-colour comparison chart providing precision in result interpretation and ease of use.

A drop of blood is placed on the strip and wiped after 30 seconds. The strip is then compared to the colour chart for a result after another 30 seconds.

BETACHEK Visual is a proven visual blood glucose test sold in over 50 countries worldwide, setting the standard for ease of use, accuracy and precision.

### BETACHEK Visual has:


- A large range from 0.5 mmol/l (9 mg/dl) to 55 mmol/l (991mg/dl).
- Excellent matching resulting from small increments in the colour chart.
- An “easy to read” pink and blue colour chart.
- Twin colour pad technology for precision in matching colours.



# BETACHEK *Visual*

## Blood Glucose Test Strips

### 3 Simple steps to monitor Blood Glucose visually



Apply blood




Wipe Blood off



Results in just 60 seconds!

Betachek v Laboratory  
Prince of Wales Hospital



Clinical study: Betachek test strips were assessed for accuracy by comparison with formal hospital laboratory method. There was good correlation between the Betachek test and the laboratory reference method ( $R=0.9777$ ).